

Antibodies Product Datasheet

Product Name:	KRR1 Antibody
Catalog Number:	C50890
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KRR1
Concentration:	1 mg/ml
Gene Name:	KRR1 HRB2
Protein Name:	KRR1 small subunit processome component homolog (HIV-1 Rev-binding protein 2) (KRR-R motif-containing protein 1) (Rev-interacting protein 1) (Rip-1)
Other Name:	
Human Swiss Prot No.:	Q13601
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KRR1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year