

## Antibodies Product Datasheet

<b>Product Name:</b>	TSSK 4 Antibody
<b>Catalog Number:</b>	C19097
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human TSSK4. AA range:211-260
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TSSK4
<b>Protein Name:</b>	Testis-specific serine/threonine-protein kinase 4
<b>Other Name:</b>	TSSK4; C14orf20; STK22E; TSSK5; Testis-specific serine/threonine-protein kinase 4; TSK-4; TSSK-4; Testis-specific kinase 4; Serine/threonine-protein kinase 22E
<b>Human Swiss Prot No.:</b>	Q6SA08
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	TSSK 4 Polyclonal Antibody detects endogenous levels of TSSK 4 protein.
<b>Storage and Stability:</b>	-20°C/1 year