
Antibodies Product Datasheet

Product Name:	β -Tubulin Monoclonal Antibody
Catalog Number:	N1379
Reactivity:	Zebrafish, Fulvidraco, Blowfish
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of β -Tubulin
Concentration:	1 mg/ml
Gene Name:	TUBB3
Protein Name:	Tubulin beta-3 chain
Other Name:	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
Human Swiss Prot No.:	Q13509
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects Zebrafish endogenous β -tubulin protein.
Storage and Stability:	-20°C/1 year