
Antibodies Product Datasheet

Product Name:	PAK1 Antibody
Catalog Number:	R12-3296
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	PAK1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PAK1.
Concentration:	1 mg/ml
Gene Name:	PAK1
Other Name:	Serine/threonine-protein kinase PAK 1; p21-activated kinase 1; PAK-1; p65-PAK; Alpha-PAK; PAK1
Human Swiss Prot No.:	Q13153
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PAK1 Antibody antibody detects endogenous levels of PAK1.
Storage and Stability:	Stable at -20°C for at least 1 year.