

Antibodies Product Datasheet

Product Name:	VPS24 Antibody
Catalog Number:	C48999
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human VPS24
Concentration:	1 mg/ml
Gene Name:	CHMP3 CGI149 NEDF VPS24 CGI-149
Protein Name:	Charged multivesicular body protein 3 (Chromatin-modifying protein 3) (Neuroendocrine differentiation factor) (Vacuolar protein sorting-associated protein 24) (hVps24)
Other Name:	
Human Swiss Prot No.:	Q9Y3E7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of VPS24 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year