

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | hnRNP Q Antibody |
| Catalog Number: | C10139 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human hnRNP Q. AA range:236-285 |
| Concentration: | 1 mg/ml |
| Gene Name: | SYNCRIP |
| Protein Name: | Heterogeneous nuclear ribonucleoprotein Q |
| Other Name: | SYNCRIP; HNRPQ; NSAP1; Heterogeneous nuclear ribonucleoprotein Q; hnRNP Q; Glycine- and tyrosine-rich RNA-binding protein; GRY-RBP; NS1-associated protein 1; Synaptotagmin-binding; cytoplasmic RNA-interacting protein |
| Human Swiss Prot No.: | O60506 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | hnRNP Q Polyclonal Antibody detects endogenous levels of hnRNP Q protein. |
| Storage and Stability: | -20°C/1 year |