

Antibodies Product Datasheet

Product Name:	FGF-8 Antibody
Catalog Number:	C30056
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FGF8. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	FGF8
Protein Name:	Fibroblast growth factor 8
Other Name:	FGF8; AIGF; Fibroblast growth factor 8; FGF-8; Androgen-induced growth factor; AIGF; Heparin-binding growth factor 8; HBGF-8
Human Swiss Prot No.:	P55075
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-8 Polyclonal Antibody detects endogenous levels of FGF-8 protein.
Storage and Stability:	-20°C/1 year