

Antibodies Product Datasheet

Product Name:	EphA4 Antibody
Catalog Number:	C40183
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EphA4
Concentration:	1 mg/ml
Gene Name:	EPHA4 HEK8 SEK TYRO1
Protein Name:	EphA4
Other Name:	Ephrin type-A receptor 4 (EC 2.7.10.1) (EPH-like kinase 8) (EK8) (hEK8) (Tyrosine-protein kinase TYRO1) (Tyrosine-protein kinase receptor SEK) P54764
Human Swiss Prot No.:	P54764
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human EphA4
Storage and Stability:	-20°C/1 year