

Antibodies Product Datasheet

Product Name:	SGLT3/SLC5A4 Antibody
Catalog Number:	C48781
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SGLT3/SLC5A4
Concentration:	1 mg/ml
Gene Name:	SLC5A4 SAAT1 SGLT2
Protein Name:	Low affinity sodium-glucose cotransporter (Sodium/glucose cotransporter 3) (Na(+)/glucose cotransporter 3) (Solute carrier family 5 member 4)
Other Name:	
Human Swiss Prot No.:	Q9NY91
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SGLT3/SLC5A4 at Human
Storage and Stability:	Store at -20° C/1 year