
Antibodies Product Datasheet

Product Name:	Tak1 Antibody
Catalog Number:	R12-2374
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Tak1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Tak1.
Concentration:	1 mg/ml
Gene Name:	MAP3K7
Other Name:	Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1
Human Swiss Prot No.:	O43318
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Tak1 antibody detects endogenous levels of Tak1.
Storage and Stability:	Stable at -20°C for at least 1 year.