

Antibodies Product Datasheet

Product Name:	SCP-3 Monoclonal Antibody
Catalog Number:	N0744
Reactivity:	Human
Applications:	WB
Dilutions:	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	
Immunogen:	Purified recombinant human SCP-3 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	SYCP3
Protein Name:	Synaptonemal complex protein 3
Other Name:	SYCP3; SCP3; Synaptonemal complex protein 3; SCP-3
Human Swiss Prot No.:	Q8IZU3
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	SCP-3 Monoclonal Antibody detects endogenous levels of SCP-3 protein.
Storage and Stability:	-20°C/1 year