
Antibodies Product Datasheet

Product Name:	PTEN mouse Monoclonal Antibody(2C10)
Catalog Number:	N1641
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of PTEN at AA range of 320-400
Concentration:	1 mg/ml
Gene Name:	PTEN
Protein Name:	PTEN
Other Name:	PTEN
Human Swiss Prot No.:	P60484
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PTEN protein detects endogenous levels of PTEN
Storage and Stability:	-20°C/1 year