

Antibodies Product Datasheet

Product Name:	CKMT2 mouse mAb
Catalog Number:	N0898
Reactivity:	Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CKMT2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ckmt2
Protein Name:	
Other Name:	CKMT 2;Basic-type mitochondrial creatine kinase;CKMT 2;CKMT2;CPK;Creatine kinase mitochondrial 2;Creatine kinase mitochondrial 2 (sarcomeric);Creatine kinase S-type; creatine kinase S-type, mitochondrial;KCRS_HUMAN;Mib CK;Mib-CK;mitochondrial; OTTHUMP0000
Human Swiss Prot No.:	P17540
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CKMT2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year