

Antibodies Product Datasheet

Product Name:	PLC γ 1 Antibody
Catalog Number:	B0079
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLCG1. AA range:736-785
Concentration:	1 mg/ml
Gene Name:	PLCG1
Protein Name:	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1
Other Name:	PLCG1; PLC1; 1-phosphatidylinositol 4; 5-bisphosphate phosphodiesterase gamma-1; PLC-148; Phosphoinositide phospholipase C-gamma-1; Phospholipase C-II; PLC-II; Phospholipase C-gamma-1; PLC-gamma-1
Human Swiss Prot No.:	P19174
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PLC γ 1 Polyclonal Antibody detects endogenous levels of PLC γ 1 protein.
Storage and Stability:	-20°C/1 year