

Antibodies Product Datasheet

Product Name:	VPS72 Antibody
Catalog Number:	C11112
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VPS72. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	VPS72
Protein Name:	Vacuolar protein sorting-associated protein 72 homolog
Other Name:	VPS72; TCFL1; YL1; Vacuolar protein sorting-associated protein 72 homolog; Protein YL-1; Transcription factor-like 1
Human Swiss Prot No.:	Q15906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	VPS72 Polyclonal Antibody detects endogenous levels of VPS72 protein.
Storage and Stability:	-20°C/1 year