

Antibodies Product Datasheet

Product Name:	DFF45 (Cleaved-Thr225) Antibody
Catalog Number:	L0483
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DFF45 (Cleaved-Thr225)
Concentration:	1 mg/ml
Gene Name:	DFFA DFF1 DFF45 H13
Protein Name:	DFF45 (Cleaved-Thr225)
Other Name:	DNA fragmentation factor subunit alpha (DNA fragmentation factor 45 kDa subunit) (DFF-45) (Inhibitor of CAD) (ICAD)
Human Swiss Prot No.:	O00273
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human DFF45 (Cleaved-Thr225, protein was cleaved amino acid sequence between 224-225)
Storage and Stability:	-20°C/1 year