

Antibodies Product Datasheet

Product Name:	CYCS Monoclonal Antibody(4B10)
Catalog Number:	N1513
Reactivity:	Human, Mouse, Rat, Chicken
Applications:	WB, IF/ICC, IHC-p
Dilutions:	WB: 1:1000-5000 IHC: 1:500-1000 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of CYCS
Concentration:	
Gene Name:	CYCS
Protein Name:	Cytochrome c
Other Name:	CYCS; CYC; Cytochrome c
Human Swiss Prot No.:	P99999
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CYCS protein.
Storage and Stability:	-20°C/1 year