

## Antibodies Product Datasheet

<b>Product Name:</b>	Cleaved-Caspase-3 p12 (D175) Antibody
<b>Catalog Number:</b>	L0105
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IF/ICC, IHC-p, ELISA
<b>Dilutions:</b>	WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Caspase 3. AA range:157-206
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CASP3
<b>Protein Name:</b>	Caspase3
<b>Other Name:</b>	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
<b>Human Swiss Prot No.:</b>	P42574
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Cleaved-Caspase-3 p12 (D175) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-3 p12 protein resulting from cleavage adjacent to D175.
<b>Storage and Stability:</b>	-20°C/1 year