Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: PKC β Antibody

Catalog Number: B8169

Reactivity: Human, Mouse, Rat

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other

applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human PKCB. AA range:622-671

Concentration: 1 mg/ml

Gene Name: PRKCB

Protein Name: Protein kinase C beta type

Other Name: PRKCB; PKCB; PKCB1; Protein kinase C beta type; PKC-B; PKC-beta

Human Swiss Prot No.: P05771

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: PKC β Polyclonal Antibody detects endogenous levels of PKC β protein.

Storage and Stability: -20°C/1 year