
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | GNL1 Antibody |
| Catalog Number: | C16070 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human GNL1. AA range:61-110 |
| Concentration: | 1 mg/ml |
| Gene Name: | GNL1 |
| Protein Name: | Guanine nucleotide-binding protein-like 1 |
| Other Name: | GNL1; HSR1; Guanine nucleotide-binding protein-like 1; GTP-binding protein HSR1 |
| Human Swiss Prot No.: | P36915 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | GNL1 Polyclonal Antibody detects endogenous levels of GNL1 protein. |
| Storage and Stability: | -20°C/1 year |