

Antibodies Product Datasheet

Product Name:	CCKBR Antibody
Catalog Number:	C40187
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CCKBR
Concentration:	1 mg/ml
Gene Name:	CCKBR
Protein Name:	Gastrin/cholecystokinin type B receptor (CCK-B receptor) (CCK-BR) (Cholecystokinin-2 receptor) (CCK2-R)
Other Name:	CCKBR; CCKRB; Gastrin/cholecystokinin type B receptor; CCK-B receptor; CCK-BR; Cholecystokinin-2 receptor; CCK2-R; YM3541
Human Swiss Prot No.:	P32239
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CCKBR protein(A214) detects endogenous levels of CCKBR
Storage and Stability:	-20°C/1 year