
Antibodies Product Datasheet

Product Name:	STAT3 Monoclonal Antibody(3F4)
Catalog Number:	N1544
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of STAT3
Concentration:	1 mg/ml
Gene Name:	STAT3
Protein Name:	Signal transducer and activator of transcription 3 (Acute-phase response factor)
Other Name:	STAT3; APRF; Signal transducer and activator of transcription 3; Acute-phase response factor
Human Swiss Prot No.:	P40763
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	STAT3 protein detects endogenous levels of STAT3
Storage and Stability:	-20°C/1 year