

Antibodies Product Datasheet

Product Name:	EIF1AX Antibody
Catalog Number:	C49350
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF1AX
Concentration:	1 mg/ml
Gene Name:	EIF1AX EIF1A EIF4C
Protein Name:	Eukaryotic translation initiation factor 1A, X-chromosomal (eIF-1A X isoform) (Eukaryotic translation initiation factor 4C) (eIF-4C)
Other Name:	
Human Swiss Prot No.:	P47813
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF1AX at Human, Mouse
Storage and Stability:	Store at -20° C/1 year