

---

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | PTH Monoclonal Antibody   |
| <b>Catalog Number:</b>        | N0540   |
| <b>Reactivity:</b>            | Human   |
| <b>Applications:</b>          | WB,ELISA  |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications. |
| <b>Modification:</b>          | Unmodified/ Total   |
| <b>Source:</b>                | Monoclonal Mouse  |
| <b>Purification:</b>          | Affinity purification   |
| <b>Immunogen:</b>             | Purified recombinant fragment of human PTH (aa1-115) expressed in E. Coli.          |
| <b>Concentration:</b>         |   |
| <b>Gene Name:</b>             | PTH   |
| <b>Protein Name:</b>          | Parathyroid hormone   |
| <b>Other Name:</b>            | PTH; Parathyroid hormone; PTH; Parathormone; Parathyrin                             |
| <b>Human Swiss Prot No.:</b>  | P01270  |
| <b>Formulation:</b>           | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.                  |
| <b>Specificity:</b>           | PTH Monoclonal Antibody detects endogenous levels of PTH protein.                   |
| <b>Storage and Stability:</b> | -20°C/1 year  |