

Antibodies Product Datasheet

Product Name:	β-Tubulin mouse mAb
Catalog Number:	N0875
Reactivity:	Human, Mouse, Monkey, Goat, Hamster, Rat, Chlamydomonas
Applications:	WB, IF/ICC
Dilutions:	wb 1:5000 icc 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Tubulin beta protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	tubb
Protein Name:	
Other Name:	beta 4 tubulin; Beta 5 tubulin; beta 1b tubulin; M40; TUBB1; TUBB5; tubulin beta 1 chain; Tubulin beta 4 chain; tubulin beta 5 chain; Tubulin beta; Tubulin beta chain; tubulin beta polypeptide.
Human Swiss Prot No.:	P07437
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Tubulin beta and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year