
Antibodies Product Datasheet

Product Name:	Tak1 (Phospho-Thr184) Antibody
Catalog Number:	P12-1057
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Tak1 (Phospho-Thr184) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Tak1 around the phosphorylation site of Thr184.
Concentration:	1 mg/ml
Gene Name:	MAP3K7
Other Name:	Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1
Human Swiss Prot No.:	O43318
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Tak1 (Phospho-Thr184) antibody detects endogenous levels of Tak1 only when phosphorylated at Thr184.
Storage and Stability:	Stable at -20°C for at least 1 year.