

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Cytokeratin 19 (ABT315) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0301
<b>Clone ID:</b>	ABT315
<b>Reactivity:</b>	Human
<b>Applications:</b>	IHC, WB
<b>Application Dilutions</b>	IHC 1:100-200, WB 1:500-1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	KRT19
<b>Gene Name:</b>	P08727
<b>Other Name:</b>	40 kDa keratin intermediate filament, CK 19, CK-19, CK19, Cytokeratin 19, Cytokeratin-19, K19, K1C19_HUMAN, K1CS, Keratin 19, Keratin type I 40 kD, Keratin type I 40kD, Keratin type I cytoskeletal 19, Keratin, type I cytoskeletal 19, Keratin, type I, 40 kd, Keratin-19, KRT19, MGC15366
<b>Isotype</b>	IgG2a, Kappa
<b>Localization:</b>	Colon/ Stomach
<b>Positive Tissue Control:</b>	Cytoplasmic, Membranous
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human CK19 protein, and shows no cross reaction with CK1, 5, 6, 7, 8, 10, 13, 14, 15, 17, 18, 20.
<b>Storage and Stability:</b>	-20°C/1 year