

Antibodies Product Datasheet

Product Name:	Anti-Human FSH (ABT172) Mouse Monoclonal Antibody
Catalog Number:	V0343
Clone ID:	ABT172
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	FSHB
Gene Name:	P01225
Other Name:	Follicle Stimulating Hormone, Follicle stimulating hormone beta polypeptide, Follicle stimulating hormone beta subunit, Follitropin beta chain, Follitropin subunit beta, FSH-B, FSH-beta, FSHB, FSHbeta
Isotype	IgG2a, Kappa
Localization:	Pituitary
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human FS H protein, and shows no cross reaction with LH, CG and TSH proteins.
Storage and Stability:	-20°C/1 year