

## Antibodies Product Datasheet

<b>Product Name:</b>	Cleaved-Cathepsin Z (L62) Antibody
<b>Catalog Number:</b>	L0189
<b>Reactivity:</b>	Human, Monkey
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human CATZ. AA range:43-92
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CTSZ
<b>Protein Name:</b>	Cathepsin Z
<b>Other Name:</b>	CTSZ; Cathepsin Z; Cathepsin P; Cathepsin X
<b>Human Swiss Prot No.:</b>	Q9UBR2
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Cleaved-Cathepsin Z (L62) Polyclonal Antibody detects endogenous levels of fragment of activated Cathepsin Z protein resulting from cleavage adjacent to L62.
<b>Storage and Stability:</b>	-20°C/1 year