

Antibodies Product Datasheet

Product Name:	TAF5L Antibody
Catalog Number:	R12-3583
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	TAF5L Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAF5L.
Concentration:	1 mg/ml
Gene Name:	TAF5L
Other Name:	TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L; PCAF-associated factor 65 beta; PAF65-beta; TAF5L; PAF65B
Human Swiss Prot No.:	O75529
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TAF5L Antibodyantibody detects endogenous levels of TAF5L.
Storage and Stability:	Stable at -20°C for at least 1 year.