

Antibodies Product Datasheet

Product Name:	BRUNOL6 Antibody
Catalog Number:	C50473
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BRUNOL6
Concentration:	1 mg/ml
Gene Name:	CELF6 BRUNOL6
Protein Name:	CUGBP Elav-like family member 6 (CELF-6) (Bruno-like protein 6) (CUG-BP- and ETR-3-like factor 6) (RNA-binding protein BRUNOL-6)
Other Name:	
Human Swiss Prot No.:	Q96J87
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRUNOL6 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year