

Antibodies Product Datasheet

Product Name:	SR-1B Antibody
Catalog Number:	C12009
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human 5-HT-1B. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	HTR1B
Protein Name:	5-hydroxytryptamine receptor 1B
Other Name:	HTR1B; HTR1DB; 5-hydroxytryptamine receptor 1B; 5-HT-1B; 5-HT1B; S12; Serotonin 1D beta receptor; 5-HT-1D-beta; Serotonin receptor 1B
Human Swiss Prot No.:	P28222
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SR-1B Polyclonal Antibody detects endogenous levels of SR-1B protein.
Storage and Stability:	-20°C/1 year