
Antibodies Product Datasheet

Product Name:	Cleaved-Factor Xa activated HC (I235) Antibody
Catalog Number:	L0126
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FA10. AA range:216-265
Concentration:	1 mg/ml
Gene Name:	F10
Protein Name:	Coagulation factor X
Other Name:	F10; Coagulation factor X; Stuart factor; Stuart-Prower factor
Human Swiss Prot No.:	P00742
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Factor Xa activated HC (I235) Polyclonal Antibody detects endogenous levels of fragment of activated Factor Xa activated HC protein resulting from cleavage adjacent to I235.
Storage and Stability:	-20°C/1 year