

Antibodies Product Datasheet

Product Name:	HNF-4 α (Acetyl-Lys106) Antibody
Catalog Number:	D12128
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Acetyl peptide from human protein at AA range: 106
Concentration:	1 mg/ml
Gene Name:	HNF4A HNF4 NR2A1 TCF14
Protein Name:	Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) (TCF-14) (Transcription factor HNF-4)
Other Name:	Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) (TCF-14) (Transcription factor HNF-4)
Human Swiss Prot No.:	P41235
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HNF-4 α when Acetyl occurs at Lys106
Storage and Stability:	-20°C/1 year