

Antibodies Product Datasheet

Product Name:	IL-4I1 Antibody
Catalog Number:	C21350
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human IL-4I1.
Concentration:	1 mg/ml
Gene Name:	IL4I1
Protein Name:	L-amino-acid oxidase
Other Name:	IL4I1; FIG1; L-amino-acid oxidase; LAAO; LAO; Interleukin-4-induced protein 1; IL4-induced protein 1; Protein Fig-1; hFIG1
Human Swiss Prot No.:	Q96RQ9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-4I1 Polyclonal Antibody detects endogenous levels of IL-4I1 protein.
Storage and Stability:	-20°C/1 year