

Antibodies Product Datasheet

Product Name:	β -actin Monoclonal Antibody(Mix)
Catalog Number:	N1403
Reactivity:	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB: 1:5000 IHC: 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of β -actin
Concentration:	1 mg/ml
Gene Name:	ACTB
Protein Name:	Actin cytoplasmic 1
Other Name:	ACTB; Actin, cytoplasmic 1; Beta-actin
Human Swiss Prot No.:	P60709
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous β -actin protein.
Storage and Stability:	-20°C/1 year