
Antibodies Product Datasheet

Product Name:	Smad4 (Phospho-Thr276) Antibody
Catalog Number:	A8420
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Smad4 (Phospho-Thr276)
Concentration:	1 mg/ml
Gene Name:	SMAD4 DPC4 MADH4
Protein Name:	Smad4 (Phospho-Thr276)
Other Name:	Mothers against decapentaplegic homolog 4 (MAD homolog 4) (Mothers against DPP homolog 4) (Deletion target in pancreatic carcinoma 4) (SMAD family member 4) (SMAD 4) (Smad4) (hSMAD4)
Human Swiss Prot No.:	Q13485
Formulation:	Liquid in PBS containing 50% glycerol, and 0.94% sodium azide.
Specificity:	This antibody detects endogenous levels of Smad4 (Phospho-Thr276) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year