

## Antibodies Product Datasheet

<b>Product Name:</b>	Toll-like Receptor 3 mouse mAb
<b>Catalog Number:</b>	N1045
<b>Reactivity:</b>	Transfected
<b>Applications:</b>	WB
<b>Dilutions:</b>	wb dilution 1:1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human Toll-Like Receptor 3 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	tlr3
<b>Protein Name:</b>	
<b>Other Name:</b>	CD283 ; CD283 antigen ; IIAE2 ; TLR 3 ; Tlr3 ; TLR3_HUMAN ; Toll Like Receptor 3 ; Toll-like receptor 3.
<b>Human Swiss Prot No.:</b>	O15455
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Transfected Only.
<b>Storage and Stability:</b>	-20°C/1 year