

Antibodies Product Datasheet

Product Name:	Smad1 Antibody
Catalog Number:	R12-3544
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	Smad1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Smad1.
Concentration:	1 mg/ml
Gene Name:	SMAD1
Other Name:	BSP-1; BSP1; Dwarf-A; Dwf-A; JV4-1; MADH1; MADR1; Mad-related protein 1; Mad1; Mothers against DPP homolog 1; Mothers against decapentaplegic homolog 1; Mothers-against-DPP-related-1; SMAD 1; Smad 1; Transforming growth factor- beta signaling protein-1
Human Swiss Prot No.:	Q15797
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Smad1 Antibody antibody detects endogenous levels of Smad1.
Storage and Stability:	Stable at -20°C for at least 1 year.