

Antibodies Product Datasheet

Product Name:	FLT3 (Phospho-Tyr589/591) Antibody
Catalog Number:	A9520
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FLT3 (Tyr589 and 591)
Concentration:	1 mg/ml
Gene Name:	FLT3 CD135 FLK2 STK1
Protein Name:	FLT3 (Tyr589/591)
Other Name:	Receptor-type tyrosine-protein kinase FLT3 (EC 2.7.10.1) (FL cytokine receptor) (Fetal liver kinase-2) (FLK-2) (Fms-like tyrosine kinase 3) (FLT-3) (Stem cell tyrosine kinase 1) (STK-1) (CD antigen CD135)
Human Swiss Prot No.:	P36888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse FLT3 (phospho-Tyr589 or 591)
Storage and Stability:	-20°C/1 year