

Antibodies Product Datasheet

Product Name:	CKI-alpha Antibody
Catalog Number:	R12-2621
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	CKI-alpha Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CKI-alpha.
Concentration:	1 mg/ml
Gene Name:	CSNK1A1
Other Name:	CK1-alpha; CKI-alpha; CSNK1A1; Casein kinase I, alpha isoform; KC1A; kinase CK1-alpha
Human Swiss Prot No.:	P48729
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CKI-alpha Antibody antibody detects endogenous levels of CKI-alpha.
Storage and Stability:	Stable at -20°C for at least 1 year.