

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | TEF-4 Antibody |
| Catalog Number: | C11087 |
| Reactivity: | Human, Mouse |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human TEAD2. AA range:71-120 |
| Concentration: | 1 mg/ml |
| Gene Name: | TEAD2 |
| Protein Name: | Transcriptional enhancer factor TEF-4 |
| Other Name: | TEAD2; TEF4; Transcriptional enhancer factor TEF-4; TEA domain family member 2; TEAD-2 |
| Human Swiss Prot No.: | Q15562 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | TEF-4 Polyclonal Antibody detects endogenous levels of TEF-4 protein. |
| Storage and Stability: | -20°C/1 year |