

Antibodies Product Datasheet

Product Name:	Bcl-10 Antibody
Catalog Number:	C13030
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BCL10. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	BCL10
Protein Name:	B-cell lymphoma/leukemia 10
Other Name:	BCL10; CIPER; CLAP; B-cell lymphoma/leukemia 10; B-cell CLL/lymphoma 10; Bcl-10; CARD-containing molecule enhancing NF-kappa-B; CARD-like apoptotic protein; hCLAP; CED-3/ICH-1 prodomain homologous E10-like regulator; CIPER; Cellular homolog
Human Swiss Prot No.:	O95999
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bcl-10 Polyclonal Antibody detects endogenous levels of Bcl-10 protein.
Storage and Stability:	-20°C/1 year