

Antibodies Product Datasheet

Product Name:	PR40A Antibody
Catalog Number:	C45645
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 450-530
Concentration:	1 mg/ml
Gene Name:	PRPF40A FBP11 FLAF1 FNBP3 HIP10 HYPA HSPC225
Protein Name:	Pre-mRNA-processing factor 40 homolog A (Fas ligand-associated factor 1) (Formin-binding protein 11) (Formin-binding protein 3) (Huntingtin yeast partner A) (Huntingtin-interacting protein 10) (HIP-10) (Huntingtin-interacting protein A) (Renal carcinoma a
Other Name:	
Human Swiss Prot No.:	O75400
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	PR40A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year