

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Oct-2 (ABT211) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0423
<b>Clone ID:</b>	ABT211
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	POU2F2
<b>Gene Name:</b>	P09086
<b>Other Name:</b>	class 2, Lymphoid restricted immunoglobulin octamer binding protein, Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2, NF A2, Oct 2, Oct-2, Octamer binding transcription factor 2, Octamer-binding protein 2, Octamer-binding transcription factor 2, OTF-2, OTF2, PO2F2_HUMAN, POU domain, POU domain class 2 transcription factor 2, POU2F2, transcription factor 2
<b>Isotype</b>	IgG2b, kappa
<b>Localization:</b>	Tonsil/ Appendix
<b>Positive Tissue Control:</b>	Nuclear
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Oct-2 protein.
<b>Storage and Stability:</b>	-20°C/1 year