

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | UQCRB Antibody |
| Catalog Number: | C48144 |
| Reactivity: | Human |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit,IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human UQCRB |
| Concentration: | 1 mg/ml |
| Gene Name: | UQCRB UQBP |
| Protein Name: | Cytochrome b-c1 complex subunit 7 (Complex III subunit 7) (Complex III subunit VII) (QP-C) (Ubiquinol-cytochrome c reductase complex 14 kDa protein) |
| Other Name: | |
| Human Swiss Prot No.: | P14927 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of UQCRB at Human |
| Storage and Stability: | Store at -20° C/1 year |