

Antibodies Product Datasheet

Product Name:	Smad2 Antibody
Catalog Number:	R12-3545
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Smad2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Smad2.
Concentration:	1 mg/ml
Gene Name:	SMAD2
Other Name:	JV18-1; MADH2; MADR2; Mad-related protein 2; Mothers against DPP homolog 2; Mothers against decapentaplegic homolog 2; Smad 2
Human Swiss Prot No.:	Q15796
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Smad2 Antibody antibody detects endogenous levels of Smad2.
Storage and Stability:	Stable at -20°C for at least 1 year.