

Antibodies Product Datasheet

Product Name:	TORC2 Phospho-Ser171 Antibody
Catalog Number:	A9815
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TORC2 Phospho-Ser171
Concentration:	1 mg/ml
Gene Name:	CRTC2 TORC2
Protein Name:	TORC2 Phospho-Ser171
Other Name:	CREB-regulated transcription coactivator 2 (Transducer of regulated cAMP response element-binding protein 2) (TORC-2) (Transducer of CREB protein 2)
Human Swiss Prot No.:	Q53ET0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.134% sodium azide.
Specificity:	This antibody detects endogenous levels of TORC2 Phospho-Ser171 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year