

Antibodies Product Datasheet

Product Name:	Casein Kinase I γ 1/2/3 Antibody
Catalog Number:	B8089
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CK-1 γ 1/2/3. AA range:229-278
Concentration:	1 mg/ml
Gene Name:	CSNK1G1/CSNK1G2/CSNK1G3
Protein Name:	Casein kinase I isoform gamma-1/2/3
Other Name:	CSNK1G1; Casein kinase I isoform gamma-1; CKI-gamma 1; CSNK1G2; CK1G2; Casein kinase I isoform gamma-2; CKI-gamma 2; CSNK1G3; Casein kinase I isoform gamma-3; CKI-gamma 3
Human Swiss Prot No.:	Q9HCP0/P78368/Q9Y6M4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Casein Kinase I γ 1/2/3 Polyclonal Antibody detects endogenous levels of Casein Kinase I γ 1/2/3 protein.
Storage and Stability:	-20°C/1 year