

Antibodies Product Datasheet

Product Name:	Bcl-6 (Acetyl-Lys379) Antibody
Catalog Number:	D12252
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Bcl-6 (Acetyl Lys379)
Concentration:	1 mg/ml
Gene Name:	BCL6 BCL5 LAZ3 ZBTB27 ZNF51
Protein Name:	Bcl-6 (Acetyl Lys379)
Other Name:	B-cell lymphoma 6 protein (BCL-6) (B-cell lymphoma 5 protein) (BCL-5) (Protein LAZ-3) (Zinc finger and BTB domain-containing protein 27) (Zinc finger protein 51)
Human Swiss Prot No.:	P41182
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Bcl-6 (Acetyl Lys379)
Storage and Stability:	-20°C/1 year